


Dr. Maisch

Any Column, Any Size, Any Media



LONGLIFE

用于制备的弹簧柱


MADE BY DR. MAISCH

CONTENT

- P 4 柱装填技术
- P 5 固定床柱填充相关的内在问题
- P 6 轴向压缩 - 使用活塞进行填充
- P 7 长寿命色谱柱硬件
- P 8 基于MODCOL弹簧柱原理的长寿命技术

- P 9 专利技术
- P 10 Reprosil手性填料
- P 11 其他制造商的填料
- P 12 3 μ m填料用于非手性SFC
- P 13 高分辨率制备色谱柱
3 μ m填料的性能与优势

- P 14 重新装填服务
- P 15 总结



**ONGLIFE
MADE BY DR. MAISCH**


来自欧洲最大的高效液相色谱柱制造商之一

柱装填技术

通过溶剂流动进行装填
轴向压缩-活塞装填

DAC – **D**ynamic **A**xial **C**ompression

SAC – **S**tatic **A**xial **C**ompression



Dilute slurry in column and packing reservoir
Bed compaction in column and packing through the column
Removal of packing reservoir; pressure released
Extrusion of stationary phase due to pressure release
Removal of extruded column bed
Closing the column (at ambient pressure)
Possibility of formation of voids under eluent flow pressure

固定床柱装填存在的内在问题


SAMPLES

- 需要消耗过多的填料
- 压力释放后填料会进一步损失
 - 装载能力较低
 - 效率较低
- 压力释放时可能床层破裂
- 在流动压力填料下，由于沿塔长度的压力下降而导致的床填料密度不均匀
- 在使用过程中形成空隙的可能导致寿命下降。

轴向压缩——活塞填装

轴向压缩克服了这些问题：

长寿命柱硬件




不需要多余的填料或导致填料损失

活塞释放压力时：填装过程中，内部填料颗粒受力均匀。


活塞压力不释放时：

- 不破坏床层
- 填充密度一致
- 不形成空隙
- 增加色谱柱的寿命

- 获得专利的高级制备硬件
- SAC和DAC模式通用
- 适用于 SFC
- 极高的性能和使用寿命





长寿命技术是基于
MODCOL弹簧柱原理




获得专利的长寿命原理

填充柱的床层长度可以通过使用half-tube柱插件来控制



可以插入垫圈弹簧单元代替静态垫片，用于动态轴向压缩机机制

- 填料储罐用于容纳稀浆。
- 填料储罐在填充完柱后移除，以尽量减少填充柱的总硬件长度。
- 活塞留在柱中。
- 压力没有释放。



专利的立柱法兰端板特殊设计，无需拆卸
活塞和释放填料压力即可关闭立柱。

Patent No: DE202018001788

DE202016000500111

CHIRAL REPROSIL MEDIA





OTHER MANUFACTURERS MEDIA

ReproSil Chiral-MIA,
5 µm, 250 mm L x 30 mm ID

ReproSil Chiral NR,
8 µm, 260 mm L x 50 mm ID

Zorbax SB-AQ,
5 µm, 70 mm L x 30 mm ID

Luna C18 (3)
10 µm, prep 250 mm L x 70 mm ID



ca. 90,000 N/m

ca. 60,000 N/m

ca. 80,000 N/m

ca. 40,000 N/m

TEST CONDITIONS

Mobile Phase: Heptan/IPA 85/15
Flow Rate: 30 ml/min
Temperature: Ambient
Pressure: 34 bar
Detector: UV @ 229 nm
Sensitivity: 0.5 mV

TEST CONDITIONS

Mobile Phase: Heptan/IPA 85/15
Flow Rate: 60 ml/min
Temperature: Ambient
Pressure: 12 bar
Detector: UV @ 254 nm
Sensitivity: 0.6 mV

TEST CONDITIONS

Mobile Phase: MeOH/H₂O 85/15
Flow Rate: 30 ml/min
Temperature: Ambient
Pressure: 24 bar
Detector: UV @ 254 nm
Sensitivity: 59.8 mV

TEST CONDITIONS

Mobile Phase: MeOH/H₂O 85/15
Flow Rate: 120 ml/min
Temperature: Ambient
Pressure: 10 bar
Detector: UV @ 254 nm
Sensitivity: 1.8 mV

Description:

Packing Material: ReproSil Chiral-MIA, 5 µm
Length: 250 mm ID: 30 mm
Shipping Solvent: Mobile Phase
Maximum Pressure: 130 bar
Hardware Type: LongLife
Frit: 2 µm
pH Range: 2.0 - 8.0

Description:

Packing Material: ReproSil Chiral-NR, 8 µm
Length: 260 mm ID: 50 mm
Shipping Solvent: Mobile Phase
Maximum Pressure: 250 bar
Hardware Type: LongLife
Frit: 2 µm
pH Range: 2.0 - 8.0

Description:

Packing Material: Zorbax SB-AQ, 5 µm
Length: 75 mm ID: 30 mm
Shipping Solvent: Mobile Phase
Maximum Pressure: 210 bar
Hardware Type: LongLife
Frit: 2 µm
pH Range: 2.0 - 8.0

Description:



Packing Material: Luna 100 C18(3), 10 µm
Length: 250 mm ID: 70 mm
Shipping Solvent: Mobile Phase
Maximum Pressure: 100 bar
Hardware Type: LongLife
Frit: 2 µm
pH Range: 2.0 - 8.0

Peak 1: Uracil
Peak 2: Toluene

3 µm MEDIA FOR ACHIRAL SFC

Reprospher 100 2-EP,
3 µm, 100 mm L x 50 mm ID

Reprospher 100 PEI,
3 µm, 100 mm L x 50 mm ID



ca. 90,000 N/m

ca. 110,000 N/m

TEST CONDITIONS

Mobile Phase: MeOH/H2O 85/15
Flow Rate: 60 ml/min
Temperature: Ambient
Pressure: 85 bar
Detector: UV @ 254 nm
Sensitivity: 21.1 mV

Description:

Packing Material: Reprospher 100 PEI 3 µm
Length: 100 mm ID: 50 mm
Shipping Solvent: Mobile Phase
Maximum Pressure: 200 bar
Hardware Type: LongLife SFC
Frit: 2 µm
pH Range: 2.0 - 8.0

TEST CONDITIONS


Mobile Phase: MeOH/H2O 85/15
Flow Rate: 60 ml/min
Temperature: Ambient
Pressure: 120 bar
Detector: UV @ 254 nm
Sensitivity: 58.7 mV

Description:

Packing Material: Reprospher 100 PEI 3 µm
Length: 100 mm ID: 50 mm
Shipping Solvent: Mobile Phase
Maximum Pressure: 200 bar
Hardware Type: LongLife SFC
Frit: 2 µm
pH Range: 2.0 - 8.0

HIGH RESOLUTION PREP CHROMATOGRAPHY
PREP COLUMN PERFORMANCE
WITH 3 µm PARTICLES

UP-SCALE



1 - Uracil 2 - Phenol 3 - N,N-Diethyl-M-Toluamide 4 - Toluene

250 x 4,6 mm

250 x 50 mm


BENEFITS OF LONGLIFE

DR. MAISCH

- Packed by piston
- Flexible bed length
- DAC and SAC mechanism
- Packing and repacking service
- Available column ID - 25, 30, 40, 50, 70
- Scalability to > 150 mm ID -Using ModCol column / Multipacker


保存每一列和重新包装

通常列硬件表示列值的重要部分：



长寿命包装的resil - pur Basic-C18 10µm ; 250 × 50毫米

节省多次重新包装的长寿命硬件。



SUMMARY

01

性能和稳定性非常高！
列尺寸更短与MoDcol相比。

02

技术
填料与MoDcol类似，但
填料后储层与色谱柱分离。

03

活塞停留在柱上

04

只能安装在 Dr.
Maisch HPLC

05

选择使用具有相同包
装的MoDcol列技术，
如果感兴趣的自包装
或直径 > 70毫米

06

LongLife在DAC和
SAC模式下可用：
25、30、40、50
和70毫米内径

Dr. Maisch

Any Column, Any Size, Any Media

Distributor:

Dr. Maisch HPLC GmbH
Beim Brückle 14
D-72119 Ammerbuch
T: +49 (0) 7073 50357
F: +49 (0) 7073 4216
www.dr-maisch.com
www.modcol.com
E-Mail: info@dr-maisch.com



PDF brochure
for download

