## 根据USP 621用ChroZen HPLC和UHPLC分析布洛 芬

- HPLC & UHPLC 应用

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### · ChroZen HPLC System

| Item            | Description                                                | Part No.   |
|-----------------|------------------------------------------------------------|------------|
| Pump            | ChroZen HPLC Quaternary Gradient Pump with Vacuum degasser | 9421011020 |
| Autosampler     | ChroZen HPLC Autosampler                                   | 5421011020 |
| Column oven     | ChroZen HPLC Column oven for Analytical scale              | 3421011020 |
| Detector        | ChroZen HPLC UV/Vis Detector with dual wavelength          | 7411011020 |
| Install. Option | HPLC Performance Kit (Without LC C18 Column)               | 1601011890 |
| CDS             | YL-Clarity software for single instrument of YCM HPLC      | 5301011000 |
|                 | Autosampler control of YL-Clarity                          | 5301011040 |
|                 | System Suitability Test of YL-Clarity                      | 5301011050 |
| Column          | C18 (4.6 mm x 250 mm, 5 $\mu$ m)                           | -          |

### · ChroZen UHPLC System

| Item            | Description                                           | Part No.   |
|-----------------|-------------------------------------------------------|------------|
| Pump            | ChroZen UHPLC Pump                                    | 9431011012 |
| Autosampler     | ChroZen UHPLC Cool/Heat Autosampler                   | 9451011013 |
| Column oven     | ChroZen UHPLC Column oven for Analytical scale        | 3531011030 |
| Detector        | ChroZen UHPLC UV/Vis Detector                         | 7521011020 |
| Install. Option | UHPLC Performance Kit (Without LC C18 Column)         | 9361011150 |
| CDS             | YL-Clarity software for single instrument of YCM HPLC | 5301011000 |
|                 | Autosampler control of YL-Clarity                     | 5301011040 |
|                 | System Suitability Test of YL-Clarity                 | 5301011050 |
| Column          | Poroshell C18 (2.1 mm x 100 mm, 2.7 $\mu$ m)          | -          |

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| 5 replicate injections | Resolution | Area RSD (%) |               |
|------------------------|------------|--------------|---------------|
|                        |            | Ibuprofen    | Valerophenone |
| USP Regulation         | $\geq 2.5$ | $\leq 2.0$   | $\leq 2.0$    |
| Result                 | 6.04       | 0.293        | 0.267         |
|                        | Pass       | Pass         | Pass          |

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| 5 replicate injections | Resolution | Area RSD (%) |               |
|------------------------|------------|--------------|---------------|
|                        |            | Ibuprofen    | Valerophenone |
| USP Regulation         | $\geq 2.5$ | $\leq 2.0$   | $\leq 2.0$    |
| Result                 | 3.42       | 0.638        | 0.582         |
|                        | Pass       | Pass         | Pass          |

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|                                     | USP <621>                                                                                                                                                                 | HPLC Condition | UHPLC Condition                     |
|-------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------|-------------------------------------|
| pH of mobile phase                  | $\pm 0.2$                                                                                                                                                                 | As specified   | As specified                        |
| Concentration of salts in buffer    | $\pm 10\%$                                                                                                                                                                | As specified   | As specified                        |
| Ratio of components in mobile phase | $\pm 30\%$ Relative<br>Cannot exceed $\pm 10\%$                                                                                                                           | As specified   | As specified                        |
| Wavelength of UV-visible detector   | $\pm 3 \text{ nm}$                                                                                                                                                        | As specified   | As specified                        |
| Column length                       | <i>Isocratic separation</i><br>; Column length(L) to particle size diameter(dp) ratio can be adjusted between -25% and 50%<br><i>Gradient separation</i><br>; Not allowed | As specified   | As specified                        |
| Column inner diameter               | Can be adjusted if the linear velocity maintained                                                                                                                         | As specified   | 2.1 mm<br>(As specified)            |
| Particle size                       | <i>Isocratic separation</i><br>; Column length(L) to particle size diameter(dp) ratio can be adjusted between -25% and 50%<br><i>Gradient separation</i><br>; Not allowed | As specified   | 2.7 $\mu\text{m}$<br>(As specified) |
| Flow rate                           | $\pm 50\%$ (Isocratic only)                                                                                                                                               | As specified   | 0.6 mL/min<br>(As specified)        |
| Injection volume                    | Can be adjusted as far as it is consistent with accepted precision, linearity, and detection limit                                                                        | As specified   | 1 $\mu\text{L}$<br>(As specified)   |
| Column temperature                  | $\pm 10^\circ\text{C}$                                                                                                                                                    | As specified   | As specified                        |

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