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
Item	Description	Part No.
Pump	ChroZen HPLC Quaternary Gradient Pump with Vacuum degasser 2 set	9421011020
Autosampler	ChroZen HPLC Autosampler	5421011020
Column Compartment	ChroZen HPLC Column Oven for Analytical scale	3421011020
Detector	ChroZen HPLC UV/Vis Detector with dual wavelength	7411011020
Install. Option	HPLC Performance Kit	1601011890
CDS	YL-Clarity software of Young In Chromass HPLC for single instrument control	5301011000
	Autosampler control of YL-Clarity	5301011040
Column	Agilent ZORBAX Eclipse XDB-C8 (4.6 mm x 150 mm, 5 μm) Agilent Eclipse Plus C18 (4.6 mm x 50 mm, 5 μm) Agilent InfinityLab Poroshell EC-C18 (2.1 mm x 250 mm, 4 μm)	-
Switching Valve	6-port 2-position valve, SS/PEEK, 40MPa , 0.4, 10-32 for ChroZen HPLC	3421012020

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C = Concentration of Vitamin B₇ in Sample (μg/mL)

a = Total Solution Volume (mL)

V_1 = Volume of Loaded Solution in Immunoaffinity Column (mL)

V_2 = Volume of Eluted Solution from Immunoaffinity Column (mL)

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- Korean standard methods for the analysis of food ingredients at trace level (Biotin)
- Simultaneous Determination of Vitamin D3 and K1 in Infant Formula by HPLC with column-switching
(Korean journal of food science and technology, Vol. 37, No; 6, 1024-1027)
- New nutritional concepts of some vitamins and minerals
(Korean Journal of Pediatrics, Vol. 48, No. 12)
- Adverse Effects on the skin by the deficiency of vitamins and minerals
(Korean Journal of Aesthetics and cosmetology, Vol. 10, No. 4)

60, Anyangcheondong-ro, Dongan-gu,
Anyang-si, Gyeonggi-do, 14042, Korea

TEL: +82-31-428-8700

FAX: +82-31-428-8787

E-mail: export@youngincm.com

Homepage: www.youngincm.com