

多环芳烃

EPA Method 610/8100

Peaks:

- | | | |
|------------------------|--------------------------|----------------------------|
| 1. Naphthalene | 8. Anthracene | 14. Benzo(k)fluoranthene |
| 2. 1-methylnaphthalene | 9. Fluoranthene | 15. Benzo(a)pyrene |
| 3. 2-methylnaphthalene | 10. Pyrene | 16. 3-methylcolantrene |
| 4. Acenaphthylene | 11. Benzo(a)anthracene | 17. Indeno(1,2,3-cd)pyrene |
| 5. Acenaphthene | 12. Chrysene | 18. Dibenzo(a,h)anthracene |
| 6. Fluorene | 13. Benzo(b)fluoranthene | 19. Benzo(ghi)perylene |
| 7. Phenanthrene | | |

Column: **MEGA-PAH - 0.25mm, 0.25 μ m, 30m**

Catalog Code: C-PAH-025-025-30

条件:

进样口: 冷柱进样

检测器: FID 315 $^{\circ}$ C

柱箱程序: 70 $^{\circ}$ C-4 $^{\circ}$ C/min-310 $^{\circ}$ C

载气: 氢气, 70kPa

