

黄油甘油三酯 (T24-T56)

Column: **MEGA-1 - 0.32mm, 0.10µm, 3m + 1m Built-In Ret. Gap**

Catalog Code: GAP-1-032-010-3-1

条件:

进样口: 冷柱进样

检测器: FID 350℃

柱箱程序: 60-340℃, 15℃/min

载气: 氢气, 30kPa

Peaks:

- | | |
|----------------|---------|
| 1. T24 | 13. T46 |
| 2. Cholesterol | 14. T48 |
| 3. T26 | 15. T50 |
| 4. T28 | 16. T52 |
| 5. T30 | 17. T54 |
| 6. T32 | 18. T56 |
| 7. T34 | |
| 8. T36 | |
| 9. T38 | |
| 10. T40 | |
| 11. T42 | |
| 12. T44 | |

T为C原子总数

