

有机氯杀虫剂

Peaks:

- | | | |
|----------------|------------------------|---------------------|
| 1. Trichlorfon | 8. Chlorpyrifos methyl | 15. Quinalphos |
| 2. Phosdrin | 9. Parathion | 16. Bromophos ethyl |
| 3. Thionazine | 10. Pirimiphos methyl | 17. Methidathion |
| 4. Ethoprophos | 11. Chlorpyrifos | 18. Ethion |
| 5. Phorate | 12. Malathion | 19. Trithion |
| 6. Fonophos | 13. Parathion | 20. Azinphos methyl |
| 7. Diazinone | 14. Pirimiphos | 21. Azinphos ethyl |

Column: **MEGA-POF 1 - 0.32mm, 0.25 μ m, 25m**

Catalog Code: C-POF1-032-025-25

条件:

进样口: 不分流240°C

检测器: FPD 140°C

柱箱程序: 100°C-5.5°C/min-230°C-30°C/min-270°C

载气: 氦气, 70kPa

