


香柠檬油

Column: **MEGA-SE54 FAST - 0.10mm, 0.10 μ m, 5m**

Catalog Code: F-SE54-010-010-5



条件:

进样口: 分流230 $^{\circ}$ C, 0.5 μ L进样量, 分流比1:250

检测器: FID 250 $^{\circ}$ C

柱温箱程序: 50 $^{\circ}$ C-250 $^{\circ}$ C, 30 $^{\circ}$ C/min

载气: 氢气, 0.5mL/min

样品制备: 溶解在环己烷中1%

峰:

1. α -蒎烯
2. β -蒎烯
3. 月桂烯
4. 柠檬烯
5. 对伞花烃
6. γ -萜品烯
7. 里哪醇
8. 乙酸芳樟酯

Acknowledgement: Prof. C.Bicchi, C.Brunelli et al. - Università di Torino, Dipartimento di Scienza e Tecnologia del Farmaco, Via Pietro Giuria, 9 - Torino - Italy.