

USP 467 溶剂残渣 (OVIs) Headspace on MEGA-624 FAST (G43 Phase)

Column: **MEGA-624 FAST - 0.10mm, 0.45 μ m, 10m**

Catalog Code: F-624-010-045-10

条件:

进样口: 分流 250 $^{\circ}$ C, 1:100分流比, 0.5mL气体进样量

检测器: FID 250 $^{\circ}$ C


柱温箱程序: 35-100 $^{\circ}$ C, 15 $^{\circ}$ C/min

载气: 氢气, 0.4mL/min

制样: 顶空制样在80 $^{\circ}$ C持续45分钟

Peaks:

1. Methylene Chloride
2. Chloroform
3. Benzene
4. 1,1,1-trichloroethylene
5. 1,4-dioxane



Analysis carried out with DANI MASTER-GC.