


Column: **MEGA-DEX DAC Beta** - 0.25mm, 0.25 μ m, 25m
Catalog Code: DEX-DAC-B-025-025-25



条件:

柱温箱程序: 50°C-225°C, 2°C/min

载气: 氢气, 80KPa, 恒压模式

进样口: 分流250°C, 分流比90mL/min, 0.5 μ L的进样量

检测器: FID 250°C

(S)-(-)- β -香茅醇 (CAS # 7540-51-4), Sigma-Aldrich
Car. # 303488

(R)-(+)- β -香茅醇 (CAS # 1117-61-9), Sigma-Aldrich
Car. # 303461

制样: 2 μ L样品溶解在1mL的正己烷中