

顶空固相微萃取鼠尾草的芳香标记样品

Column: MEGA-WAX 20M MS - 0.20mm, 0.20 μ m, 50m

Catalog Code: MS-WAX20M-020-020-50


条件:

进样口: 分流 230 $^{\circ}$ C, 1:20分流比

检测器: MSD EI模式(70eV),扫描范围35-350amu, 延压40ms,
离子源温度230 $^{\circ}$ C, 传输管线温度280 $^{\circ}$ C

柱温箱程序: 50 $^{\circ}$ C, 5 $^{\circ}$ C/min

载气: 氦气, 1mL/min



Acknowledgments: many thanks to Barbara Sgorbini, Cecilia Cagliero, Chiara Cordero and Carlo Bicchi, Università di Torino, Dipartimento di Scienza e Tecnologia del Farmaco, Via P. Giuria, 9, Torino, Italy.

Full reference paper: Journal Of Chromatography A, 1376 (2015) 9-17.